


# Lbs of LN2 to go From Elevated Temperature to Ambient



Ambient  
Temp test  
day  
22 degree  
(C)

Temp.

Temp (Scale Left)

Lbs LN2 (Scale Right)

Temperature Degrees (C)

LN2

LN2

LN2

Temp.

LN  
2

LN2  
Temperature


0:00:00	0:01:39	0:02:36	0:03:18	0:03:48	0:04:15	0:04:42	0:05:03	0:05:24	0:05:45	0:06:06
1.7	1.46	1.28	1.02	0.88	0.74	0.54	0.48	0.34	0.2	0
118.876	115.718	105.421	95.245	85.7	75.941	65.117	55.676	45.761	35.647	25.116

# Lbs of LN2 to go From Ambient to Lowered Temperature

Ambient  
Temp Test  
day  
22 degree  
(C)

Temp (Scale left)

Lbs LN2 (Scale Right)



LN2	0	0.16	0.18	0.32	0.46	0.52	0.82	1.04	1.14	1.42	1.96	1.98	2.02
Temperature	25.116	15.458	5.348	-5.501	-15.163	-24.764	-35.591	-44.98	-54.689	-65.63	-65.835	-65.894	-65.957